



# Cell type: Primary Chondrocytes

Performance of Eppendorf Cell Culture Consumables (TC treated)

**Procedure:** Human primary articular chondrocytes were seeded in appropriate culture medium (CGM™ BulletKit™). Cell Culture vessels were incubated at 37 °C with a humidified 5 % CO<sub>2</sub> atmosphere. Cell growth and morphology was monitored on day 2, 3 and 4 post-seeding.



Bars represent average count of two wells of an Eppendorf 6-Well Cell Culture Plate



Representative light microscope images (100X magnification)

The tissue culture (TC) treated surface of Eppendorf Cell Culture Consumables allows optimal attachment and homogeneous growth of **primary chondrocytes**.

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf AG · 22331 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)