

# Eppendorf Certificate

## Certificate of Quality

### Eppendorf Serological Pipets

The Eppendorf Quality Management System is certified according to ISO 9001, ISO 13485, and ISO 14001. This results for our customers to a constantly high and reliable quality of products and processes covering the complete product life time. Products are manufactured in a class 8 clean room which meets the ISO 14644-1 specifications.

The raw material, that is used to produce the Eppendorf Serological Pipets, complies with USP Class VI standards. Eppendorf confirms that Eppendorf Serological Pipets and the packaging materials comply to REACH regulation of the European Union EC No 1907/2006 amended by Commission Regulation EC No. 552/2009.

Material suppliers do not use or intentionally incorporate the following agents into the materials used for the production of Eppendorf Serological Pipets:

- Slip agents (including oleamide, erucamide, stearamides)
- Biocides (including di(2-hydroxyethyl)methyl dodecyl ammonium salts (DiHEMDA))
- Plasticizers (including phthalates)
- Bisphenol A
- Latex

Eppendorf Serological Pipets have been tested for in-vitro cytotoxicity according to ISO 10993-5 and passed its requirements.

For certification of absence of detectable DNA, DNase, and RNase, and pyrogens as well as for sterilization, please refer to lot-specific certificates of the respective product. Lot-specific purity certificates can be downloaded on [www.eppendorf.com/certificates](http://www.eppendorf.com/certificates).

January, 2020



Joana Tziolis  
Expert Liquid Handling  
LH-Consumables

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)



Monika Schneider  
Head of Global Quality Management &  
Regulatory Affairs

ISO 9001  
Certified

ISO 13485  
Certified

ISO 14001  
Certified